



DIASAFE[®]*plus* 

... offering the **plus** in value

Ultrapure Dialysis Fluid - Good Dialysis Practice



DIASAFE®*plus* Fluid Filter



DIASAFE®*plus* Fluid Filter

- New filter design concept
 - improved flow distribution/header design
 - hygienic connectors: DIAFIX® Lock System, exclusive to DIASAFE®*plus*
- fast and easy installation/exchange
- Increased surface area - 2.2 m²
 - increased safety
 - increased retention capacity
- Filter life-time*:
 - 100 treatments or 12 weeks
 - Improved chemical resistance



* integrity test passed



Special Filter Holder: Handling

The DIAFIX™ Lock System facilitates the simple, rapid and safe connection of the filters to the machine:




1. Open locking
2. Insert DIASAFE *plus*
3. Close locking

3 rapid handling steps only




DIASAFE[®]plus



Dialysis fluid quality

- Amount of dialysis fluid contacting the patients blood: 360 L / week
(HD, QD = 500 ml/min, 3 x 4h)



- HDF: Higher QD
Substitution fluid (directly injected into extracorporeal blood circuit)
- Removal of inorganic contaminations → water treatment system
Microbiological quality → Ultrafilters



Backtransport in dialysers

Backdiffusion + Backfiltration = Backtransport

LPS

Lipid A fragments
other LPS fragments


2000 – 4000 Da
< 8000 Da

Peptidoglycan
Muramin peptides

400 – 1000 Da

β_2 - M

11 800 D




Membranes function by size exclusion principles only.

→ *The dialysis membrane is no “one-way street”*




Backtransport in dialysers



Clinical advantages of using ultrapure dialysis fluid

- Reduction of chronic micro-inflammation oxidative stress
- Lower inflammatory response (cytokine expression, c-reactive protein synthesis)
- Less hypoalbuminemia, as a consequence of lower inflammatory response
- Improved iron utilization and better response to EPO
→ better control of anaemia
- Reduction of plasma pentosidine levels (AGE compound)
- Higher quality treatment short- and long-term
- Lower incidence of dialysis-related Amyloidosis



Dialysis Fluid

- Components (water, concentrates) according to chemical and microbiological standards.
- Microbiological quality defined for dialysis fluid before filtration to increase safety:

	Water* European Pharmacopoeia 2005	Concentrate EN 13867 (2002) Ph. Eur. 2005	Dialysis Fluid	ONLINE substitution fluid (ISO Standard (draft))**
Microbial counts (CFU/mL)	<100	<100	<100	<10 ⁻⁶
Endotoxin levels (IU/mL)	<0.25	<0.5	<0.25	<0.01 (below detection limit)

* Ph. Eur. 5.0: haemodialysis solutions, concentrated, water for diluting (reverse osmosis)

** ISO Standard for dialysis fluid for haemodialysis and related therapies (draft)

- Maintenance and disinfection of all systems dialysis fluid comes into contact with is necessary!



Ultrapure Fluid

Incoming dialysis fluid filtered by two DIASAFE[®]*plus*- filters
→ resulting in ultrapure dialysis fluid:

Effective retention of microorganisms and endotoxins by DIASAFE[®]*plus* (Weber et al. 2000)

		Bacteria count [CFU/mL]	Endotoxins [IU/mL]
Effective retention by 1 st DIASAFE [®] <i>plus</i>	Dialysis fluid	< 0.01	< 0.01
Effective retention by 2 nd DIASAFE [®] <i>plus</i>	Dialysis fluid	0	< 0.01

→ Filtration across DIASAFE[®]*plus* effectively removes bacteria and endotoxins



Features

- high retention capacity for bacterial products
- cross-flow filtration
- hygienic connector technique
- new filter design
- filter life-time
- chemical resistance



Benefits

- ultrapure dialysis fluid
- no accumulation of contaminants
- no source of contamination
- high filtration rates at low pressure
- 100 treatments or 12 weeks
- disinfection and cleaning agents validated for HD machine





for filtration of dialysis fluids
and for
preparation of substitute solution
for ONLINE-HDF and HF treatments



Specifications DIASAFE[®] *plus*

Filtration rate	5 mL/min mmHg (3.75 L/min bar)
max. filtration pressure	2 bar
Membrane material	Fresenius Polysulfone [®]
Membrane surface	2.2 m ²
Housing material	Polypropylene
Potting material	Polyurethane
Sterilisation	INLINE Steam
Article No.	500 8201

The logo features a large, stylized, handwritten-style script of the word "Polysulfone" in black. Below the script, the word "Fresenius" is written in a smaller, blue, sans-serif font on the left, and "Polysulfone[®]" is written in a smaller, blue, sans-serif font on the right.



Specification Diasafe[®] *plus*

**Installation/
Filter exchange**

**max. 100 Treatments,
resp. 12 Weeks**

**Disinfection /
Cleaning**


**Citrosteril, 84 °C
Puristeril
Diasteril**

Sporotal, 11x

Treatment and cleaning cycles with Sporotal are automatically counted.




Competitor Comparison: Dialysate Fluid Filter System



The ULTRA system by Gambro, with its set-up of ultrafilters, provides ultrapure dialysis fluid and non-pyrogenic substitution fluid for HDF or HF.

Gambro: **3 Ultrafilters** for preparation of substitution fluid for HDF or HF



↔ FMC: only 2 DIASAFEplus filters necessary



Competitor Comparison: Dialysate Fluid Filter System

Machine	Dialysate Fluid Filter	Membrane	Disinfectants	Duration of Use
FMC 4008 / 5008	2 DIASAFE <i>plus</i>	Fresenius Polysulfon	Puristeril Citrosteril Diasteril	100 treatments resp. 3 months
Gambro	2 U8000S + 1 U2000S for ONLINE- HDF	Polyethersulfon	Peracetic acid Citric acid	4 weeks



Summary

DIASAFE® *plus*...

...ensures the safe production of ultrapure dialysis and substitution fluid.

...has excellent endotoxin-retention capabilities of its Fresenius Polysulfone® fibres.


And in addition:

- DIASAFE® *plus* can be used for a maximum of 12 weeks.
- The DIAFIX™ Lock System facilitates simple, rapid and safe connection of the filters to the machine.




New Brochure

- Clinical advantages of using ultrapure dialysis fluid
- Cardioprotective Haemodialysis
→ ultrapure dialysis fluid supports reduction of chronic inflammation
- Ultrapure dialysis fluid is a mandatory prerequisite for ONLINE HDF
→ major selling point of the 5008



New Literature Compendium

Art no. 733 971 1



Thank you for your attention !

